

INCREASE INSTRUMENT RELIABILITY



Buy 2 Get 1 Free

2009 Fall Special

EST Analytical
Preventative Maintenance Services

Preventative maintenance can help reduce instrument downtime and increase reliability. Instrument failure wastes time, sample, effort and money. A regularly scheduled program of preventative maintenance is important to your success.

EST Analytical recognizes the value of preventative maintenance and would like to provide these services at a discount for a limited time. Through the month of October, purchase two PM Services and receive the third at no charge. (Equal or lesser value)



MSD	Agilent 597X	\$ 975
GC	Agilent 5890/6890/7890	\$ 685
P&T Concentrator	EST Analytical / OI / Tekmar	\$ 500
P&T Autosampler	EST Analytical / Archon	\$ 1545

Pricing includes Parts and Labor

MSD PM Includes

- Source cleaning
- Vacuum maintenance
- Clean electronics bay
- Pump down / leak check
- Review of MSD condition
- Applications / Method Review

GC PM Includes

- Inlet maintenance
- Column maintenance
- Clean electronics bay
- Check oven and flapper motor
- Chemical test
- Applications / Method Review

Concentrator PM Includes

- Trap replacement
- Water Mngt replacement (Where applicable)
- Valve testing
- Clean electronics bay and fans
- Application / Method Review

Autosampler PM Includes

- XYZ Rebuild
- Bushings and bearings
- Gripper maintenance
- Syringe maintenance
- Qualification / Inspection
- Application / Method Review

Schedule your service today by calling:

513.642.0111

503 Commercial Drive
Fairfield, OH 45014



EST
analytical