

[View this message in a browser.](#)



New Applications  
Get 10% off Parts  
Orders!  
New Air Analysis  
Systems



[Forward this email](#)

## New Application Note - VPH Analysis



The analysis of Volatile Petroleum Hydrocarbons (VPH) requires purge and trap concentration in tandem with a Gas Chromatograph (GC) for compound separation. For compound determination, a Flame Ionization Detector and a Photoionization Detector (FID/PID) are used in series. Finally, a Mass Spectrometer (MS) is also an option for confirmatory analysis. The method is employed for the analysis of ground water, soils and/or any other solid waste that can be extracted in methanol. To download this application, [click here](#). To see a complete list of application notes, visit our [application section](#) of our website.

---

## New Online Store - Now save 10%!



Save 10% when you order \$300 or more now until May 31st! Just enter coupon code **EST010** when checking out to get your discount!

EST Analytical is your source for all your purge and



trap parts. We strive to make doing business with EST enjoyable and have created a new online store to make buying even easier. To visit our store and create an account, go to [www.estanalyticalstore.com](http://www.estanalyticalstore.com).

We carry an extensive inventory for all brands, and our tech support team will help you identify the

parts you need.

Let us know what you think!

---

## New Cannister System for TO-14 and TO-15 Analysis



If you are doing TO 14, TO 15 analysis you need to talk to us!

After years of listening to requests from customers for EST to enter this market, EST now has a complete canister system. Backed by our team of 15 field service engineers an application support headed up by a PhD chemist with more than 25 years or air analysis experience, we

want to be your solution for air analysis. For more information, [visit us](#).



---

## Get 20% Off Your Next Service Call!

Use our service team to service your GC, GC/MS, or purge and trap products in the month of April or May and get 20% off on labor and travel for your second service call (parts not included). For more information or to set up a service call, contact our [Technical Support Team](#).



## New Online Forum



In order to increase productivity and help customers troubleshoot systems, we have created a new [online forum](#) section on our web site. Post questions or review answers to previously posted issues. Can't find what you are looking for, try our new [online chat](#) to talk directly with our world class tech support team.

---

Visit us on Facebook



[www.estanalytical.com](http://www.estanalytical.com)

[Unsubscribe or update your email address.](#)



